



NORTH
AMERICAN
DRÄGER

[RETURN TO TECHNICAL SERVICE BULLETIN TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Technical Service Department
3122 Commerce Drive
Telford, Pennsylvania 18969
Telephone: (215) 721-5402
(800) 543-5047
Facsimile: (215) 723-5935

TECHNICAL SERVICE BULLETIN #152

To: All Qualified Service Organizations

From: Ron Greco

Date: December 30, 1996

Product: All Narkomedes with one piece table tops.

Change: The table top has been modified with a relief "slot" to the inside corner where it meets the flowmeter shield.

Rationale: Flowmeter shields were cracking at contact point.

Timing: The next Service or PMS.

Procedure: Inspect the inside corner of the table top where it meets the flowmeter shield. If it appears the shield is touching the table top, proceed with the following steps using the attached drawing as a reference.

Mark the flowmeter shield with a removable marker where the table top and the flowmeter shield meet. Remove the table top from the machine. With a fine tooth file, make a chamfer from the mark down to the bottom of the shield. This distance is approximately 1/4 inch. This process should be repeated until the flowmeter shield does not touch the table top at contact point.

Cost: There is no cost to the customer.

Parts: None

If you have any questions, please contact Technical Support by phone at 800-543-5047 or by fax at 215-723-5935.

